


Razek Arthroscopy Cannulas

INDICATIONS

Razek Arthroscopy Cannulas are devices designed for the creation and maintenance of access routes in the joints. They allow the insertion or removal of instruments, implants and sutures in arthroscopic surgical procedures.

n.° ANVISA: 80356130017
Merely illustrative images






TECHNICAL SPECIFICATIONS

- Sterilization: ethylene oxide;
- Sterilization validity: 2 years;
- Single-use product.


MODELS



Cannula 6,5 mm x 70 mm (740040000)
Threaded cannula with plug measuring 6,5 x 70 mm.



Cannula 8,5 mm x 70 mm (740120000)
Threaded cannula with plug measuring 8,5 x 70 mm.



Cannula 8,5 mm x 90 mm (740050000)
Threaded cannula with plug measuring 8,5 x 90 mm.



www.razek.com.br
+55 16 2107.2345