






# Maxcut

**PERCUTANEOUS METHOD  
FOR SPINE SURGERY**





**DIFFERENT TYPES OF ACTIVE TIPS**



## INDICATIONS

The Maxcut kit provides all the disposables needed to perform disc decompression. Its Bone Debridement Cannula removes structures, enabling access to the disc. The Puncture Cannula performs the initial puncture and through it the Nitinol Wire is properly positioned to later conduct the non-disposable instruments, always by the same path. The 330 mm Ablator is responsible for dissection, hemostasis and fulguration in places of difficult access, thanks to its actuation system and its flexible tip that provides greater precision in movements. The Irrigation Equipment promotes continuous fluid flow, minimizing possible thermal damage to adjacent tissues as well as associated bleeding.

Merely illustrative images

## MODELS

### The Maxcut kit consists of the following items:

- 01 unit of 330 mm Ablator
- 01 unit of Irrigation Equipment  
(Perforator - Evolution - STK - ATX)
- 01 unit of Puncture Cannula
- 01 unit of Nitinol Wire
- 01 unit of Bone Debridement Cannula  
Types of active tips: OBC - TBPC - TBC - DBC - RBC - AMCC - INCC  
Types of fittings: R - A - F - L - E - M - W - K



www.razek.com.br  
+55 16 2107.2345