


Razek Hip Ablator

INDICATIONS

The Razek Hip Ablator is used in monopolar arthroscopic surgical procedures. It is meant for ablation and coagulation of tissues in the hip.

Its flexible tip provides directional control of the end of your stem, making it easy to access hard-to-reach regions.


Merely illustrative images



Flexible Tip that provides directional control of the end

TECHNICAL SPECIFICATIONS

- Sterilization: ethylene oxide;
- Validity of sterilization: 2 years;
- Single-use product.



MODELS

200 Razek Hip Ablator - 880470000

Stem 200 mm long

210 Razek Hip Ablator - 880480000

Stem 210 mm long

220 Razek Hip Ablator - 880490000

Stem 220 mm long

230 Razek Hip Ablator - 880510000

Stem 230 mm long

240 Razek Hip Ablator - 880520000

Stem 240 mm long

250 Razek Hip Ablator - 880530000

Stem 250 mm long

330 Razek Hip Ablator - 881390000

Stem 330 mm long



www.razek.com.br

+55 16 2107.2345