


Razek Arthroscopy Cannulas

INDICATIONS

Razek Arthroscopy Cannulas are devices designed for the creation and maintenance of portals in the joints. They allow the insertion or removal of instruments, implants and sutures in arthroscopic surgical procedures.

Merely illustrative images





TECHNICAL SPECIFICATIONS

- Sterilization: ethylene oxide;
- Sterilization validity: 2 years;
- Single-use product.

MODELS



Cannula 6.5 mm x 70 mm (740040000)
Threaded cannula with obturator measuring 6.5 x 70 mm.



Cannula 8.5 mm x 70 mm (740120000)
Threaded cannula with obturator measuring 8.5 x 70 mm.



Cannula 8.5 mm x 90 mm (740050000)
Threaded cannula with obturator measuring 8.5 x 90 mm.



Cannula 8.5 mm x 110 mm (741940000)
Threaded cannula with obturator measuring 8.5 x 110 mm.



www.razek.com.br



+55 16 2107 2345



+55 16 2107 2346

