

Suture Compression Hook

INDICATIONS

The Suture Compression Hook is a disposable device designed to allow manipulation of the suture through the soft tissue in arthroscopic surgical procedures.





FitPass (742640000)



SCP (741960000)



FS (742620000)



EXS2 (741980000)



LIN (742000000)



SCM (741990000)



SCM2 (742660000)



BIP (742010000)



SCPIII (742610000)



STS I (742670000)



STS II (742680000)



EXS1 (741970000)




FP (742420000)



FP Hook +
FP Hook Inserter

FitPass

Razek Arthroscopy Instrument



Merely illustrative images

The Razek Arthroscopy Instrument - **FitPass** provides accurate suture passage in rotator cuff repair. Easy suture capture and retrieval with a single movement.

FitPass suture passer is safe and easy to use.



CHARACTERISTICS

- Suture Passer tip with a width of 4.8 mm, compatible with cannulas with diameters of 7.5 mm or larger;
- Lengths: 170 mm;
- Code: 742464100.



Giving value to what really matters



www.razek.com.br



+55 16 2107 2345



+55 16 2107 2346

