

KneePass

The Razek Arthroscopy Instrument - **KneePass** provides precise suture passage in meniscus repair, with easy capture and retrieval of the wire in a single motion.


KneePass Suture Passer is safe and easy to use.




See more details!

Merely illustrative images

**USE SIMULATION
Vertical Meniscus
Repair**



KneePass - R086



LOW PROFILE


Its geometry provides a safe passage between the meniscus and the femoral condyle, allowing the suturing of different injuries.

TECHNICAL FEATURES

- Suture Passer tip with a width of 4.5 mm;
- Length: 137mm;
- Code: 742464400.

RZK

Suture Passer Hook



TECHNICAL FEATURES

- Used in conjunction with the KneePass when passing sutures;
- Has Nitinol Passer;
- Code: 742700000.



Giving value to what really matters



www.razek.com.br



+55 16 2107 2345



+55 16 2107 2346

