


STP Pain Plus Cannula Set


INDICATIONS

The STP Pain Plus Cannula Set is a single-use set of instruments that allows for safe and accurate peripheral nerve stimulation and anesthetic or drug block in spine, shoulder, hip, and knee procedures. It has an insulating Nano-coating over its entire length, leaving only the bevel exposed, ensuring a precise stimulus in the desired region. Its high-intensity echogenic tip allows precise localization in ultrasound-guided procedures.


Merely illustrative images



STP PAIN CANNULA



EXTENDER TUBE



MODELS

▪ **STP Pain Plus Cannula Set 50 (Code: 921210100)**

This kit comes with 2 units of STP Pain Cannula 50 and 1 unit of Extender Tube

▪ **STP Pain Plus Cannula Set 80 (Code: 921210200)**

This kit comes with 2 units of STP Pain Cannula 80 and 1 unit of Extender Tube

▪ **STP Pain Plus Cannula Set 100 (Code: 921210300)**

This kit comes with 2 units of STP Pain Cannula 100 and 1 unit of Extender Tube

▪ **STP Pain Plus Cannula Set 120 (Code: 921210400)**

This kit comes with 2 units of STP Pain Cannula 120 and 1 unit of Extender Tube

STP Pain Cannula



Merely illustrative images

MODELS

▪ STP Pain Cannula 50 (Code: 921200100)

This model features 1 cannula with 50 mm length



▪ STP Pain Cannula 80 (Code: 921200200)

This model features 1 cannula with 80 mm length



▪ STP Pain Cannula 100 (Code: 921200300)

This model features 1 cannula with 100 mm length



▪ STP Pain Cannula 120 (Code: 921200400)

This model features 1 cannula with 120 mm length

