

Razek Aspiration Tip


The Disposable Tip manufactured by Razek is used in monopolar surgical procedures. It is indicated for tissue ablation, coagulation, and cauterization. It features an integrated aspiration system.

The Razek Aspiration Tip eliminates the need for an electro-surgical pen, as it is an integral part of the product.

Merely illustrative images

**SURGICAL PEN WITH
ERGONOMIC
HANDLE**

**With an
integrated
system for
suction of
bubbles from
the process**





Ergonomic Handle

+ comfort and safety during use

COAG activation button

CUT activation button

Polymer insulated

TECHNICAL SPECIFICATIONS

- Sterilization: ethylene oxide
- Sterilization Validity: 2 years
- Single-use product

TWO DIFFERENT TYPES OF ACTIVE TIPS



45°
Code: 881520000



90°
Code: 881520100

Active tip with ceramic coating

